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1	INTRODUCTION	20
1.1	General Information	20
1.2	Symbols used	20
2	INTENDED USE	21
2.1	Scope of Application	21
3	PRODUCT DESCRIPTION.....	21
3.1	General Information.....	21
3.2	Holding System (Product).....	22
3.3	Head Components.....	22
3.4	Articulated Arms	22
3.5	Columns.....	23
3.6	Bases (Product).....	23
4	INSTALLATION AND USE	23
4.1	Receiving Inspection	23
4.2	Warranty.....	23
4.3	Assembly	24
4.4	Disassembly	24
4.5	Operation Instructions for Head Components	25
4.6	Operation Instructions for Articulated Arms.....	26
4.7	Operation Instructions for Bases	26
5	REPROCESSING.....	26
5.1	General Requirements.....	26
5.2	Pre-Treatment	27
5.3	Machine Cleaning and Disinfection	28
5.4	Sterilization	29
5.5	Maintenance	30
5.6	Storage	30
5.7	Re-Usability.....	30
5.8	Obligation to Notify	31
5.9	Disposal	31

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1 Introduction

Thank you for buying our product. You have chosen a high quality Swiss product.

1.1 General Information

Please read these instructions carefully and keep them in a safe place. Improper handling can cause harm to the patient or damage the products. Use universal precautions when handling contaminated or biohazardous components / materials.

Our delivery and mailing address is:

Baitella AG
Thurgauerstrasse 70
CH-8050 Zurich
Switzerland

Phone: +41 44 305 80 00
Fax: +41 44 305 80 05
E-Mail: info@baitella.com
Internet: www.fisso.com

Copyrights:

All rights reserved. Duplication, adaptation, or translation, of any part of this document, without previous written authorization by Baitella AG is prohibited, except within the framework of the copyright regulations.

Reserved rights to technical changes!

The illustrations and specifications in this instruction for use may differ slightly from actual products.

1.2 Symbols used

Symbol	Definition	Danger	Consequence
	DANGER!	Immediate danger for people	Death or serious injuries
	WARNING!	Possible danger for people or objects	Health damages or severe material damages
	ATTENTION!	Possible danger to objects	Material damages

Additional Symbols used:

Symbol	Description
	Symbol for products that are put into use according to the regulation (EU) 2017/745 on medical devices
	Steam or dry-heat sterilization
	Note Additional hints or important information
	See instruction for use
	Manufacturer
	Distributor
	Medical device (in this instruction for use indicated as «product»)
	Article number

	Date of manufacture
LOT	Lot number
SN	Serial number
QTY	Quantity
Rx only	Usage only by trained medical personnel

2 Intended Use

2.1 Scope of Application

This product is used to hold and position instruments, endoscopes / optics and other devices. It must only be used by trained medical personnel capable of judging and controlling any danger to patients. If not the case, the user assumes all responsibility.

DANGER!

Unauthorized changes or modifications to the FISSO holding system are prohibited for safety reasons.

DANGER!

All FISSO components are designed for optimal performance when used together with the FISSO holding system. If a product / component from another manufacturer is used, the user takes full responsibility.

DANGER!

The product will be delivered **non sterile**. Prior to the initial and each following use, the product has to be cleaned, disinfected and sterilized as well as checked for visible irregularities and malfunctioning according to the indications given in this manual.

Guideline for medical devices (conformity declaration)

This product corresponds to the regulation (EU) 2017/745 on medical devices.

For a „**non-grounded**“ patient positioning (complying to the safety norms EN IEC 60601-1 and 2-2, VDE 0750 Part 2 and VDE 0753) an insulated rail clamping base must be used (Art. No 5.240 with a separate instruction for use).

Special care should be taken when these products are used together with **high-frequency applications**. Contacts between the product and the high-frequency device must be avoided.

3 Product Description


3.1 General Information

The following sections will describe the products / components based on examples. The product range along with all related information can be reviewed on www.fisso.com.

3.2 Holding System (Product)

A complete holding system consists of:

- head component
- articulated arm (double-arm system: two articulated arms)
- column
- base

All surgical products / components can be steam-sterilized up to 134 °C.

3.3 Head Components

Head components vary depending on their suitability for the following applications:

- **Combined Instrument holder**
e.g. order no. 5.422

Is intended for the secure, continuously variable clamping of instruments.

- **Varioholder**
e.g. order no. 5.716

Is intended for the secure, continuously variable clamping of round instruments.

- **Optic Holder**
e.g. order no. 5.447

Serves to securely clamp instruments using a form-fitting construction that does not harm or damage endoscopes / optics


Amendment for Sales in the United States:

US Federal Law restricts the endoscope holder system to sale by or on the order of a physician. Indications for Use: FISSO Holding System consists of a table-mounted endoscope holder system intended for use by surgeons to hold endoscopes and arthroscopes with a diameter of 4mm to 10mm during general diagnostic and therapeutic procedures. The device is also intended for use by qualified surgeons for holding endoscopes during diagnostic and therapeutic neurologic procedures.

In addition, head component designs vary depending on the following methods of fixation:

- **Permanently attached to the Articulated Arm**

- **Interchangeable** (QR = Quick release system)
A quick release system consists of: a quick release holder (Pos. 1, permanently attached to the articulated arm) and a quick release adapter (Pos. 2, permanently attached to the head component).

3.4 Articulated Arms

The FISSO articulated arms can be positioned exactly in any desired position. All three joints can be fixed and quickly released with the central handle (Z). The purely mechanical central clamping unit allows easy and quick positioning of endoscopes / optics and other instruments.

3.5 Columns

The columns serve to adjust the height of the articulated arm. The columns are permanently attached to the articulated arms and the column end takes the form of either a conical end piece or a stop lid.

Columns vary depending on their design:

- **Column, straight**
e.g. order no. 5.462

- **Column, L-Shape**
e.g. order no. 5.495

- **Column, U-Shape**
e.g. order no. 5.496

Column end:


- **End Piece**
e.g. order no. 93.685

- **Stop Lid**
e.g. order no. 93.689


3.6 Bases (Product)

The bases clamp the columns and are affixed to the corresponding location, e.g. on the rail of an operating table (rail clamping base). The bases vary depending on their suitability for the following applications:


- **Rail Clamping base**
e.g. order no. 5.280

Serves to clamp onto rails. The column height is adjustable.

- **Radial setting Clamp**
e.g. order no. 5.750

Serves to clamp onto rails. The column height and tilt angle are adjustable.

- **Universal Base**
e.g. order no. 5.251

Serves to clamp onto rails, round tubes or plates (depending on the design). The column height is adjustable.

For a „**non-grounded**“ patient positioning (complying to the safety norms EN IEC 60601-1 and 2-2, VDE 0750 Part 2 and VDE 0753) an insulated rail clamping base must be used (Art. No 5.240 with a separate user manual).

4 Installation and Use

4.1 Receiving Inspection

Check the product immediately after receipt for eventual transport damages and completeness. Complaints can only be considered if the seller or freight forwarder is immediately notified. In this case, there must immediately send a damage protocol to the next FISSO representative or to the company Baitella AG.

4.2 Warranty

The manufacturer assumes warranty for materials and manufacturing defects as well as proper functioning of a product for the duration of 24 months starting with the delivery of products to purchaser. Not covered by the warranty are normal wear and the consequences of inappropriate handling, combinations with other products, or damage caused by the user. Other or further claims, in particular any claims for damage are excluded.

The warranty applies only if the Baitella AG repair service is used. Please also note the General Sales and Delivery Conditions.

4.3 Assembly

**WARNING!**

All handling must be done with as little force as possible and as much force as necessary!

**WARNING!**

The assembly and handling of the products must be done manually and without additional tools.

1. Attach the base to the rail of the surgical table and secure it by turning the handle clockwise from below (Pos. 1, see figure section 3.6). Pay attention during attachment that the upper hook (Pos. 3, see figure section 3.6) grabs behind the rail.
2. If necessary remove the stop lid (see section 3.5) from the column by turning counter clockwise. Place the column with one hand into the column opening of the base and secure the column with the other hand by turning the handle (Pos. 2, see figure section 3.6) clockwise. Then screw the stop lid into the column clockwise and tighten.
3. The desired height and orientation of the column including the articulated arm can be adjusted at any time by loosening the handle (Pos. 2, see figure section 3.6).
Important: Hold the column with one hand and loosen the handle with the other hand. When the desired position is reached, tighten the handle.

**DANGER!**

If the base or the column is not correctly fixed, these components may become loose and may cause injuries to the patient.

**ATTENTION!**

The handle (Pos. 2, see figure section 3.6) should not be tightened at any time without the column first being in place.


4.4 Disassembly

The disassembly of the holding set must be performed also without additional tools and occurs in the reverse order of the assembly:


1. If necessary remove the stop lid (see section 3.5) from the column by turning counter clockwise.
Important: Hold the articulated arm with one hand and loosen the handle (Pos. 2, see figure section 3.6) with the other hand by turning the handle counter clockwise.
2. Pull the articulated arm together with the column out of the opening of the base. Thread the stop lid lightly into the column to prevent loss.
3. Remove the base from the rail of the surgical table after loosening the handle (Pos. 1, see figure section 3.6) by turning the handle counter clockwise.

4.5 Operation Instructions for Head Components


- **Combined Instrument Holder** e.g. order no. 5.422

- **Varioholder** e.g. order no. 5.716

- **Optic Holder** e.g. order no. 5.447

By turning the handle (Pos. 1) counter clockwise the instrument holder is opened and the instrument or wound hook can be placed between the clamping jaws.
By turning the handle clockwise the device is secured.

By turning the handle (Pos. 1) clockwise (in the direction of  symbol) the varioholder is opened and the round instrument can be placed between the clamping jaws. By turning the handle (Pos. 1) counter clockwise (in the direction of  symbol), the continuously variable bracket (Pos. 2) closes and the round instrument is secured in place.

By turning the handle (Pos. 1) counter clockwise the optic holder is opened and the endoscope / optic can be inserted from the front.
By turning the handle clockwise the endoscope / optic is secured

WARNING (Optic Holder)!


The endoscope / optic can be damaged if the handle is tightened too strongly.

WARNING (Optic Holder)!

Do not try to change the position of the endoscope / optic by turning the camera. This may damage the endoscope / optic. Always open the central handle (Z) of the articulated arm for repositioning (see section 4.6)

- **Quick Release System (QR)**

Pushing the button (Pos.1) on the quick release holder disengages the locking mechanism and the head component can either be inserted for assembly or be removed again for disassembly. Hold the button (Pos. 1) with one hand and operate the head component with the other hand.

DANGER!

Ensure the proper assembly of the quick release adapter (Pos. 2) by listening for the audible latching of the locking mechanism.

DANGER!

Unintentionally pressing the button (Pos. 1) disengages the quick release adapter (Pos. 2) and along with it the head component. The button (Pos. 1) must not be pushed in whenever the system is in use.


Amendment for sales in the United States:

US Federal Law restricts the endoscope holder system to sale by or on the order of a physician.

Indications for Use: FISSO Holding System consists of a table-mounted endoscope holder system intended for use by surgeons to hold endoscopes and arthroscopes with a diameter of 4mm to 10mm during general diagnostic and therapeutic procedures. The device is also intended for use by qualified surgeons for holding endoscopes during diagnostic and therapeutic neurologic procedures.

4.6 Operation Instructions for Articulated Arms

Holding-Releasing-Positioning-Fixing-Controlling:

1. Hold the articulated arm with one hand in the anterior segment and manipulate the central handle (Z) with the other hand.
2. To release turn the central handle (Z) counter clockwise as far as necessary.
3. Move the articulated arm into the desired position.
4. To secure, turn the central handle (Z) clockwise.
5. Controlling: Check that the articulated arm is firmly tightened and functioning properly.

DANGER!

- If the articulated arm is not tightened properly it can become loose and move, potentially causing injuries.
- Do not put too much weight on the articulated arm.
- The articulated arm can transmit electrical current and heat to the patient. Avoid the contact between the articulated arm and any sources of electrical current or heat.

DANGER!

- A damaged product can lead to serious injuries.
- Only use products in perfect condition and check functionality.

ATTENTION!

- The fixation of the articulated arm is based on the principle of friction. Changing the position without loosening the clamping mechanism can cause damage and will shorten the lifespan of the articulated arm.
- The articulated arm can be adjusted with little force. If the central handle (Z) has been fully loosened it has to be turned in the clockwise direction!

ATTENTION!

Residues of physiological saline solution (e.g. Sodium Chloride) affect the metallic surface. Saline solution must not come into contact with the product nor enter the articulated arm.

4.7 Operation Instructions for Bases

The general instructions for use for the base are described in section 4.3 Assembly and section 4.4 Disassembly.

DANGER!

FISSO bases are designed for optimal performance when used together with the FISSO holding system. If a base from another manufacturer is used, the user takes full responsibility.

- Radial setting clamp order no. 5.750

WARNING!

Close the gear coupling in such a manner, that the teeth grasp into each other and that they do not stand on each other

5 Reprocessing

5.1 General Requirements

The product must be cleaned, disinfected, and sterilized before each use; this specifically applies also for the first use after delivery because all products are delivered **non-sterile** (cleaning and disinfection after removing the transport protection package; sterilization after packaging). Effective cleaning and disinfection is an absolute pre-condition for effective sterilization.

The user is responsible for the sterility of the product and that only sufficiently validated instrument- and product-specific procedures are used for cleaning, disinfection, or sterilization, that the apparatus used (disinfection apparatus, sterilization apparatus) are serviced and checked regularly, and that the validated parameters are respected in each cycle.

Please consider additionally the laws as well as the hygiene regulations of the medical practice or hospital. This applies specifically also to the differing requirements regarding effective prion inactivation (not applicable for USA).

**NOTE**

All products will be delivered non-sterile. Please clean, disinfect, and sterilize the products before the first use according to the instructions.

**WARNING!**

All transport packaging and other packaging, protective caps and foils, etc. must be removed before using the products and should be disposed ecologically.

5.2 Pre-Treatment

Directly after use (maximum 2 hours) major debris must be removed from the products.

Tighten the central handle (Z) of the articulated arm and place it under running water or a disinfection liquid. The disinfecting agent should be aldehyde-free (otherwise blood soiling will set), should have proven efficacy (e.g. VAH/DGHM- or FDA-licensed, or CE-marking), should be suitable to disinfect the product and compatible with the product (see section material durability in chapter 5.3).

Please note that the disinfectant used for the pre-treatment serves only to protect people and cannot replace the mandatory disinfection step that has to be performed subsequently after cleaning.

**WARNING!**

With the exception of the central handle (Z) on the articulated arm, all handles of the products must be in the open position during pre-treatment.

**WARNING!**

The central handle (Z) of the articulated arm must be tightened during pre-treatment so that impurities cannot enter the arm.

**WARNING!**

If the head component is equipped with a quick release system (QR, see section 3.3), the head component (QR) must be removed during pre-treatment.

**WARNING!**

The handle and the fixation screw of the radial setting clamp (Pos. 1 and 2, see figure section 3.6) must be removed during pre-treatment.

**WARNING!**

The articulated arm must not be immersed in liquid.

**ATTENTION!**

Damage due to inappropriate cleaning! Use as much cleaning agent and cleaning liquid as necessary and remove surplus cleaning agent with a dry cloth.

**ATTENTION!**

Components made of aluminium alloys will be damaged by basic (pH > 9) cleaning agents and solvents.

**ATTENTION!**

Residues of physiological saline solution (e.g. Sodium Chloride) affect the metallic surfaces. Remove residues of physiological saline solution with a cloth wetted in clear water and dry the product with a lint-free cloth.

**ATTENTION!**

Only use a soft brush or a clean soft cloth that are used only for this purpose for manual removal of debris, never use metal brushes or steel wool.

**NOTE**

Always observe the instruction for use from the manufacturer of the cleaning and disinfection solutions.

5.3 Machine Cleaning and Disinfection

Cleaning and disinfection of the products must be performed in a disinfecting apparatus. A manual procedure must not be used.

**WARNING!**

With the exception of the central handle (Z) on the articulated arm, all handles of the products must be in the open position during cleaning / disinfection.

**WARNING!**

The central handle (Z) of the articulated arm must be tightened during machine cleaning / disinfection.

**WARNING!**

If the head component is equipped with a quick release system (QR, see section 3.3), the head component (QR) must be removed during machine cleaning / disinfection.

**WARNING!**

The handle and the fixation screw of the radial setting clamp (Pos. 1 and 2, see figure section 3.6) must be removed during machine cleaning/disinfection.

When choosing a cleaning and disinfecting apparatus make sure that:

- the cleaning and disinfecting apparatus has a proven efficacy (e.g. DGHM- or FDA-approval/clearance, or CE-marking according to DIN EN ISO 15883),
- a proven program for thermal disinfection (A0-value > 3000 or – in case of an older apparatus – minimally 5 min at 90°C) is used whenever possible (when using chemical disinfection there is a risk of traces of disinfection agent remaining on the product),
- the program used is suitable for the product and contains sufficient rinsing cycles,
- only sterile or germ-poor (maximum 10 germs/ml) water is used for final rinse,
- the air used for drying is filtered and
- the cleaning and disinfecting apparatus is serviced and checked regularly.

When choosing a cleaning agent, make sure that:

- it is basically suitable for cleaning the product,
- if no thermal disinfection is used, appropriate disinfecting agent with proven efficacy (e.g. VAH/DGHM- or FDA-approval/clearance, or CE-marking) is additionally used and that this is compatible with the cleaning agent used and
- the chemicals used are compatible with the products (see also Material Durability).

**NOTE**

The instructions for use of the cleaning and disinfecting agent manufacturer must be followed.

Procedure:

1. Place the products with tightened central handle (Z) in the cleaning and disinfecting apparatus. Make sure that the products do not touch each other.
2. Start the program.
3. Remove the products from the cleaning and disinfecting apparatus after the program has ended.
4. Control and package the products as soon as possible. If necessary let them dry first on a clean surface (see following section control and sterilization packaging).

Proof for the general suitability of this procedure for an effective machine cleaning and disinfection was provided by an independent accredited test laboratory using the cleaning and disinfecting apparatus G7836CD (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh, Germany) and the cleaning agent Neodisher medizym (Dr. Weigert GmbH & Co. KG, Hamburg, Germany). The procedure described above was applied.

Material Durability

When choosing cleaning and disinfecting agents, please make sure that they do not contain the following:

- organic/inorganic, or oxidizing acids (pH should not be less than 4, the use of neutral/enzymatic agents is recommended),
- bases (pH should not be higher than 9, the use of neutral/enzymatic agents is recommended),
- Solvents (alcohol, benzene, acetone etc.),
- Phenol,
- Chlorine, bromine, iodine,
- Chlorine salts (in particular ammonium chloride compounds), chlorinated/halogenated hydrocarbons and
- Oxidizing agents, peroxide, hypochlorite.

Control

Check after cleaning, respectively cleaning/disinfection the functionality of all components of the product to ensure smooth operation throughout the intended range of motion, corrosion, damaged surfaces, chipping, and dirt and sort out damaged products (limitation of re-use, see section 5.7 Re-usability). Products that are still dirty must be cleaned and disinfected again.

5.4 Sterilization

Sterilization Packaging

Package the product in single-use sterilization packaging material (single or double-packaging material) and/or sterilization containers that meet the following requirements:

- according to DIN EN ISO/ANSI AAMI ISO11607 (additionally DIN EN 868-8 Sterilization container),
- suitable for steam-sterilization (temperature resistance up to 138°C (280°F), sufficient permeability to steam), for USA: With FDA clearance,
- sufficient protection of the instruments and of the sterilization packaging against mechanical damage and
- regular checked sterilization container.

Steam Sterilization

Only the following sterilization procedures must be used:

- Fractionated vacuum procedure¹ (with sufficient drying of the product²)
- Steam-sterilization apparatus according to DIN EN 13060, DIN EN 285 and AAMI ST 79 (for USA: With FDA clearance) respectively
- Validated according to DIN EN ISO/ANSI AAMI ISO 17665 (valid commissioning (IQ/OQ) and product-specific performance evaluation (PQ))
- Maximum sterilization temperature 134°C (273°F); plus tolerance according to DIN EN ISO/ANSI AAMI ISO 17665
- Sterilization time (time of exposure to the sterilizing temperature) minimally 20 min. at 121°C (250°F), or alternatively minimally 4 min³ at 132 / 134°C (270 / 273°F), but radial setting clamp (order no. 5.750) must only be sterilized minimally 4 min³ at 132 / 134°C (270 / 273°F)

¹ The less efficient gravitation procedure may be used only if the fractionated vacuum procedure is not available. Additional product-, sterilization apparatus-, and procedure-specific validation must be performed under the responsibility of the user (if necessary considerably longer sterilization times should be applied).

² The required drying time depends directly on parameters which are in the sole responsibility of the user (configuration and density of the loading, condition of the steam-sterilization apparatus ...) and must be determined therefore by the user. Nevertheless drying times below 20 min. should be avoided.

³ 18 min. for prion inactivation

**WARNING!**

All handles of the products and the central handle (Z) of the articulated arm must be in the open position during sterilization.

**WARNING!**

By releasing the central handle (Z) the articulated arm may suddenly slide downwards and can cause injury or damage. To avoid this put the articulated arm down, hold it with one hand and with the other hand release the central handle (Z).

**WARNING!**

If the head component is equipped with a quick release system (QR, see section 3.3), the head component must be removed during sterilization.

**WARNING!**

The handle and the fixation screw of the radial setting clamp (Pos. 1 and 2, see figure section 3.6) must be inserted loose during sterilization.

**ATTENTION!**

The products must not be exposed at temperatures higher than 138°C (280°F)!

Proof for the general suitability of this procedure to effectively sterilize the products was provided by an independent accredited test laboratory using the steam sterilization apparatus EuroSelectomat (MMM Münchener Medizin Mechanik GmbH, Planegg, Germany) respectively HST 6x6x6 (Zirbus technology GmbH, Bad Grund, Germany) and by using a fractionated vacuum procedure. The procedure described above was applied.

Never use the flash-sterilization procedure. Additionally do not use hot air sterilization, irradiation sterilization, formaldehyde or ethylene oxide sterilization, or plasma sterilization.

5.5 Maintenance

Special maintenance is not required. The product should be sent back to the manufacturer or distributor if the action of moving parts or the clamping force is impaired, or in case of damage.

Instrument oil or grease may never be used.

**NOTE**

Repair must be carried out by Baitella AG. If not respected, the warranty becomes void!

**ATTENTION!**

Please send back only reprocessed products to the distributor / manufacture for repair (cleaned, disinfected, and sterilized).

5.6 Storage

**WARNING!**

Store the articulated arms such that they are not damaged.

- Sterilized products must be completely dry prior to storage.
- Sterile, packed products should be stored in a designated, limited access area that is well ventilated, and provides protection from dust, moisture, insects, vermin and temperature and humidity extremes.
- Sterile packaged products should be examined carefully prior to opening to ensure that there has been
- no loss of package integrity.

5.7 Re-Usability

The products can be re-used up to 500 cycles if they are not damaged and have been cleaned and sterilized according to the instructions provided. Every re-use beyond the 500 cycles and the use of a damaged or dirty product are in the responsibility of the user. If a damaged product is used or if a product is re-used without having been cleaned and sterilized, then the company is not liable for any damage causes.

In case of disrespect all liability is excluded.

5.8 Obligation to Notify

The purchaser shall inform the manufacturer about risks and recalls, as well as about reportable incidents, insofar as a product supplied by manufacturer is affected.

5.9 Disposal

The products must be correctly disposed of in accordance with national regulations and medical guidelines.

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Alle rechten voorbehouden. Duplicatie, aanpassing of vertaling van enig deel van dit document, zonder voorafgaande schriftelijke toestemming van Baitella AG is verboden, behalve in het kader van het copyright.
Reparationer måste utföras av Baitella AG. Om detta inte respekteras kommer garantin att upphöra att gälla!
Novid Surgical, LLC