

Positioning Examples



NOVID
SURGICAL

**FISSO Patient
Arm Positioner**



NOVID
SURGICAL

27 Siemon Company Drive
Suite 104W
Watertown, CT 06795
Telephone 860.274.5900
www.novidsurgical.com
info@novidsurgical.com


Representative



FISSO
Swiss Made

Quick, reliable set-up

Our innovative FISSO patient arm positioner is a surgical table accessory with a versatility that has no equal. Providing the user with reliable, one step ease of use means less time in positioning and more time for more important clinical tasks. The heart of our system is our signature FISSO arm assembly which features three rotation joints; when the desired position is reached the entire assembly can be locked in place from the central knob.



Viscoelastic pad

- Ideal anatomical form
- Sealed construction for optimal cleaning
- Highest patient comfort and safety

Large spherical joints

- Smooth rotation and positioning capabilities
- Exceptional holding strength
- Robust design ensures superior performance

Central locking knob

- Ergonomic design for ease of use
- Maximum leverage for tightening and loosening
- Smooth rounded surfaces



6646.01
Lateral
Arm Support



6656.02
Vertical
Arm Support