


Simply the best fixation.



FISSO
Swiss Made

Quick release system for various surgical applications



Features:

- Interchangeable quick release head components allow various surgical applications (MIS, open surgery etc.)
- Maximum Flexibility: Just change the quick release head component according to your application
- Easy, fast and secure exchange of quick release head component with pushbutton
- Functional design of the quick release provides worry free mounting and ergonomic handling
- Quick release holder interface guarantees full compatibility for present and future quick release head component applications
- Steam sterilizable (134° C)

Simply the best fixation.



FISO
Swiss Made

Quick release system Surgery

Set includes two quick release head components 5.716 and 5.731



3840.77

Order no.	Description	Head component Order no.	Articulated arm Order no. Length	Column Order no. Length	Base Order no. Type
3840.77	Two quick release head components (QR), quick release holder (QR), L-column, radial setting clamp	5.716 5.731	2.760 400 mm	5.495-05 450×100 mm	5.750 Radial setting clamp


QR = Quick Release
Subject to modifications

Quick release head component

Designed for various surgical applications



5.716



5.731

Order no.	Description	Application	Connection	Range	Weight kg
5.716	Varioholder (QR)	Round instruments	Quick release adapter	Ø 4.5–12.5 mm	0.18
5.731	Combined instrument holder (QR)	Instruments	Quick release adapter	Ø 4–18 mm Opening 0–18 mm	0.19

QR = Quick Release
Subject to modifications

Novid Surgical, LLC
27 Siemon Company Drive, Suite 104W
Watertown, CT 06795
Phone (860) 274-5900, Fax (860) 274-5906
info@novidsurgical.com, www.novidsurgical.com



BMPKe 03/17