


## Articulated Arms

We are happy to assist you in determining the ideal articulated arm for your specific application. Below are examples from our extensive product portfolio.

### Articulated arms Non-sterilizable

#### Bottom thread

For fixing columns  
and bases



#### Top thread

For fixing holders  
and adapters


Article #	Type	Length mm	Bottom thread	Top thr	Knob	Weight kg
7.330-VE	3600-28	280	M12x1	M8x1.25	plastic, grey	0.78
7.335-VE	3600-28 S	280	M12x1	M8x1.25	plastic, grey	0.78
7.340-VE	3600-40	400	M12x1	M8x1.25	plastic, grey	0.9
7.345-VE	3600-40 S	400	M12x1	M8x1.25	plastic, grey	0.86
7.365-VE	3600-55 S	550	M12x1	M8x1.25	plastic, grey	0.97
2.510-VE	3600-70	700	M12x1	M8x1.25	plastic, grey	1.25
2.610-VE	3600-90	900	M12x1	M8x1.25	plastic, grey	1.40
3.910-VE	3500-90	900	M12x1	-	plastic, grey	1.03
7.440-VE	4600-40	400	M12x1	M8x1.25	plastic, grey	1.90
7.460-VE	4600-55	550	M12x1	M8x1.25	plastic, grey	2.05
7.470-VE	4600-70	730	M12x1	M8x1.25	plastic, grey	2.29
2.212-VE	4600-100	1000	M12x1	M8x1.25	plastic, grey	2.70
7.530-VE	5600-32	320	M12x1	M12x1	Aluminium, red	2.50
7.550-VE	5600-46	460	M12x1	M12x1	Aluminium, red	2.70
7.570-VE	5600-70	730	M12x1	M12x1	Aluminium, red	3.30
2.260-VE	6600-32	320	M12x1	M12x1	Aluminium, red	2.93
2.270-VE	6600-46	460	M12x1	M12x1	Aluminium, red	3.13
2.275-VE	6600-56	560	M12x1	M12x1	Aluminium, red	3.33
2.280-VE	6600-70	730	M12x1	M12x1	Aluminium, red	3.58

## Articulated Arms

We are happy to assist you in determining the ideal articulated arm for your specific application. Below are examples from our extensive product portfolio.

### Articulated arms

#### 134°C Sterilizable



Article #	Type	Length mm	Bottom thread	Top thread	Knob	Weight kg
4.051-VE	XS-13	130	M6x1	M6x1	Stainless steel	0.25
4.060-VE	XS-19	190	M6x1	M6x1	Stainless steel	0.30
2.860-VE	3800-28	280	M12x1	M8x1.25	Aluminium, red	0.80
2.760-VE	3800-40	400	M12x1	M8x1.25	Aluminium, red	0.90
2.790-VE	3800-55	550	M12x1	M8x1.25	Aluminium, red	1.00
3.530-VE	4800-40	400	M12x1	M8x1.25	Aluminium, red	1.85
3.520-VE	4800-55	550	M12x1	M8x1.25	Aluminium, red	2.00
2.101-VE	5800-50	560	M12x1.75	M10x1.5	Plastic, black	3.45