


Articulated Arms

We are happy to assist you in determining the ideal articulated arm for your specific application. Below are examples from our extensive product portfolio.

Articulated arms

134°C Sterilizable



Article #	Type	Length mm	Bottom thread	Top thread	Knob	Weight kg
4.051-VE	XS-13	130	M6x1	M6x1	Stainless steel	0.25
4.060-VE	XS-19	190	M6x1	M6x1	Stainless steel	0.30
2.860-VE	3800-28	280	M12x1	M8x1.25	Aluminium, red	0.80
2.760-VE	3800-40	400	M12x1	M8x1.25	Aluminium, red	0.90
2.790-VE	3800-55	550	M12x1	M8x1.25	Aluminium, red	1.00
3.530-VE	4800-40	400	M12x1	M8x1.25	Aluminium, red	1.85
3.520-VE	4800-55	550	M12x1	M8x1.25	Aluminium, red	2.00
2.101-VE	5800-50	560	M12x1.75	M10x1.5	Plastic, black	3.45


Articulated Arms

We are happy to assist you in determining the ideal articulated arm for your specific application. Below are examples from our extensive product portfolio.

Articulated arms Non-sterilizable

Bottom thread

For fixing columns
and bases



Top thread

For fixing holders
and adapters

Article #	Type	Length mm	Bottom thread	Top thr	Knob	Weight kg
7.330-VE	3600-28	280	M12x1	M8x1.25	plastic, grey	0.78
7.335-VE	3600-28 S	280	M12x1	M8x1.25	plastic, grey	0.78
7.340-VE	3600-40	400	M12x1	M8x1.25	plastic, grey	0.9
7.345-VE	3600-40 S	400	M12x1	M8x1.25	plastic, grey	0.86
7.365-VE	3600-55 S	550	M12x1	M8x1.25	plastic, grey	0.97
2.510-VE	3600-70	700	M12x1	M8x1.25	plastic, grey	1.25
2.610-VE	3600-90	900	M12x1	M8x1.25	plastic, grey	1.40
3.910-VE	3500-90	900	M12x1	-	plastic, grey	1.03
7.440-VE	4600-40	400	M12x1	M8x1.25	plastic, grey	1.90
7.460-VE	4600-55	550	M12x1	M8x1.25	plastic, grey	2.05
7.470-VE	4600-70	730	M12x1	M8x1.25	plastic, grey	2.29
2.212-VE	4600-100	1000	M12x1	M8x1.25	plastic, grey	2.70
7.530-VE	5600-32	320	M12x1	M12x1	Aluminium, red	2.50
7.550-VE	5600-46	460	M12x1	M12x1	Aluminium, red	2.70
7.570-VE	5600-70	730	M12x1	M12x1	Aluminium, red	3.30
2.260-VE	6600-32	320	M12x1	M12x1	Aluminium, red	2.93
2.270-VE	6600-46	460	M12x1	M12x1	Aluminium, red	3.13
2.275-VE	6600-56	560	M12x1	M12x1	Aluminium, red	3.33
2.280-VE	6600-70	730	M12x1	M12x1	Aluminium, red	3.58