


Patient positioning / arm support system



Features:

« one step set up »

Central clamping mechanism for easy, fast and secure positioning of the arm support

Rapid locking base allows connection to the side rails of the OP table

Uncompromising, efficient handling and safe holding

Space saving on the OR table, compact design

Anatomically shaped visco-elastic (VE) – pad

Simply the best fixation.



FISO
Swiss Made

Patient arm support lateral



6646.01

Order no.	Description	Pad	Articulated arm Order no. Length	Carrying capacity	Base Order no.	Weight
6646.01	Arm support lateral 6600-46, strong model	visco-elastic	2.270 460 mm	10 kg	5.751	4.30 kg

Subject to modifications

Patient arm support vertical



6656.02

Order no.	Description	Pad	Articulated arm Order no. Length	Carrying capacity	Base Order no.	Weight
6656.02	Arm support vertical 6600-56, strong model	visco-elastic	2.275 560 mm	6 kg	5.750	5.40 kg

Subject to modifications

Additional components



93.283

Order no.	Description	Length	Width	Quantity	Weight
93.283	Fastening straps	600 mm	50 mm	2 pcs	0.06 kg

Subject to modifications



27 Siemon Company Drive Suite 104W
Watertown, CT 06795 Telephone 860.274.5900
www.novidsurgical.com info@novidsurgical.com



BMPase 02/18