

## FISSO Anesthesia

The FISSO holding systems with the central rapid fixation was developed in close cooperation with anesthesiologists. The purely mechanical central clamping unit enables quick and space-saving positioning of breathing tubes. Safe positioning is guaranteed throughout the anesthesia.



Art. # 3628.13

### Anaesthetic tube holder

circuit holder plate, articulated arm 280 mm, column, base

#### Consisting of:

- circuit holder plate for various  $\varnothing$ : 14.5/22/23.5 mm, coaxial >30 mm (Art.# 5.407)
- articulated arm 280 mm (Art.# 7.330)
- column 50 mm (Art.# 5.480)
- rails clamp Fisso (Art.# 5.260-02)



Art. # 3640.13

### Anaesthetic tube holder

circuit holder plate, articulated arm 400 mm, column, base

#### Consisting of:

- Circuit holder plate for various  $\varnothing$ : 14.5/22/23.5 mm, coaxial >30 mm (Art.# 5.407)
- Articulated arm 400 mm (Art.# 7.340)
- Column 50 mm (Art.# 5.480)
- Rails clamp Fisso (Art.# 5.260-02)



Art. # 2628.10

### Anaesthetic tube holder light

circuit holder plate, articulated arm 280 mm, column, base

#### Consisting of:

- Circuit holder plate for various  $\varnothing$ : 14.5/22/23.5 mm, coaxial >30 mm (Art.# 5.407)
- Articulated arm 280 mm (Art.# 3.800)
- Column 50 mm (Art.# 5.479)
- Rails clamp Fisso (Art.# 5.261)