


Simply the best fixation.

Holder with double arm

without rails clamp and quick release head component (QR)



Order no.	Description	Head component Clamping range	Articulated arm Length	Column Length	Weight
3840.48	Instrument holder double arm, L-column	5.422 (2x) Ø 4–18 mm, 0–18 mm	2.760 (2x) 400 mm	82.660 500x242 mm	4.10 kg
3840.98	Quick release holder double arm (QR), L-column	5.700 (2x)	2.760 (2x) 400 mm	82.660 500x242 mm	4.10 kg

Subject to modifications, QR=Quick Release

Quick release head component (QR)

Designed for various surgical applications




Order no.	Description	Application	Connection	Range	Weight
5.716	Varioholder (QR)	Round instruments	Quick release adapter	Ø 4.5–12.5 mm	0.18 kg
5.731	Combined instrument holder (QR)	Instruments	Quick release adapter	Ø 4–18 mm Opening 0–18 mm	0.19 kg

Subject to modifications, QR=Quick Release

Rails clamp

Designed for OR table standard rails



Order no.	Description	Material	Rail dimension	Column	Weight
5.752-VE	Radial rails clamp	Stainless steel	Height: 25–32 mm Width: 9–10 mm	For revolving and height-adjustable column Ø 16 mm	0.90 kg
5.780-VE	Rails clamp	Stainless steel	Height: 25–32 mm Width: 9–10 mm	For height-adjustable column Ø 16 mm	0.90 kg

Subject to modifications