




FISSO
Swiss Made

Simply the best fixation

Instrument holder without rails clamp



Order no.	Description	Head component Clamping range	Articulated arm Length	Column Length	Weight
3840.40	Instrument holder, straight column	5.422 Ø 4-18 mm, 0-18 mm	2.760 400 mm	5.468-05 400 mm	1.70 kg
3840.41	Instrument holder, L-column	5.422 Ø 4-18 mm, 0-18 mm	2.760 400 mm	5.495-05 450x100 mm	2.00 kg
3840.42	Instrument holder, U-column	5.422 Ø 4-18 mm, 0-18 mm	2.760 400 mm	5.496 266x330x116 mm	2.15 kg
3840.55	Instrument holder, straight column	5.715 Ø 4.5 mm-12.5 mm	2.760 400 mm	5.468-05 400 mm	1.70 kg

Subject to modifications

Rails clamp



5.752-VE



5.780-VE

Order no.	Description	Material	Rail dimension	Column	Weight
5.752-VE	Radial rails clamp	Stainless steel	Height: 25–32 mm Width: 9–10 mm	For revolving and height-adjustable column Ø 16 mm	0.90 kg
5.780-VE	Rails clamp	Stainless steel	Height: 25–32 mm Width: 9–10 mm	For height-adjustable column Ø 16 mm	0.90 kg

Subject to modifications



27 Siemon Company Drive
Suite 104W
Watertown, CT 06795 Telephone 860.274.5900
www.novidsurgical.com info@novidsurgical.com



BMPSe 10/20