

Universal Radial Clamp for US, UK and European Surgical Table Rails



Easy to mount. Easy to use. Exceptional strength.



The Universal Radial Clamp

- Infinite rotation
- Fits all table rails
- Accepts 16 - 20mm accessory posts
- Steam sterilizable
- Height adjustable

Order #10-308-16 Radial Clamp
for US, UK and European rails
16mm (5/8") hole for accessory post.

Order #10-308 Radial Clamp
for US, UK and European rails
20mm hole for posts 16-20mm



Novid Surgical is the proud US Distributor for Reison AB line of surgical table accessories.

Novid Surgical, LLC

27 Siemon Company Drive, Suite 104W, Watertown, CT 06795

Tel: 860-274-5900, Fax: 860-274-5906 www.novidsurgical.com