



Simply the best fixation.

Ventilation



3670.40

3590.30

Order no.	Description	Head component	Articulated arm Length	Column Length	Rail Clamp	Rail dimension	Weight
3670.40	Respiratory Circuit Holder	Twin 5.428-02 Single 5.427-02	2.510 700 mm	5.462 300 mm	5.260-02	Height: 25-35 mm Width: 8-10 mm	2.10 kg
3590.30	Respiratory Circuit Holder Light	Twin 5.428-02 Single 5.427-02	3.910 900 mm	5.462 300 mm	5.261-02	Height: 25-35 mm Width: 8-10 mm	1.85 kg

Subject to modifications

Head components



5.428-02-VE



5.427-02-VE



5.790-VE



5.433-02-VE

Order no.	Description	Range	Material	Weight
5.428-02-VE	Twin Tube Holder	for double tubes \varnothing 15-27 mm	PA	0.02 kg
5.427-02-VE	Single Tube Holder	for coaxial tubes \varnothing 24-33 mm	PA	0.02 kg
5.790-VE	Twin Tube Holder	for double tubes \varnothing 15 / 21.5 / 28 mm	Alu	0.70 kg
5.433-02-VE	Filter / Tube Support	for \varnothing 22 / 27 mm, length 250 mm	PA / Alu	0.08 kg

Subject to modifications

Anaesthesia



3628.13

3640.13

Order no.	Description	Head component	Articulated arm Length	Rail Clamp	Rail dimension	Weight
3628.13	Anaesthetic Tube Holder	tube plate for \varnothing 14.5 / 22 / 23.5 mm	7.330 280 mm	5.260-02-VE	Height: 25-35 mm Width: 8-10 mm	1.22 kg
3640.13	Anaesthetic Tube Holder	tube plate for \varnothing 14.5 / 22 / 23.5 mm	7.340 400 mm	5.260-02-VE	Height: 25-35 mm Width: 8-10 mm	1.33 kg

Subject to modifications

