

Patient Positioning / Arm Support System



Supine Position



Beach Chair Position



Side Position



Prone Position

Features:

- « one step set up »
- Central clamping mechanism for easy, fast and secure positioning of the arm support
- Rapid locking base allows connection to the side rails of the OP table
- Uncompromising, efficient handling and safe holding
- Space saving on the OR table, compact design
- Anatomically shaped visco-elastic (VE) – pad



Patient Arm Support Lateral



6646.01

Order no.	Description	Pad Order no.	Articulated arm Order no.	Carrying capacity	Rail Clamp Order no.	Rail dimension	Weight
6646.03	Arm support lateral 6600-46, strong model incl. Fastening Straps 5.890-VE (2 pcs.)	visco-elastic 93.279-VE	Length 2.270-VE 460 mm	10 kg	5.751-VE	Height: 25-35 mm Width: 9-10 mm	4.30 kg

Subject to modifications

Patient Arm Support Vertical



6656.02

Order no.
6656.04

	Description	Pad Order no.	Articulated arm Order no.	Carrying capacity	Rail Clamp Order no.	Rail dimension	Weight
Subject to modifications	Arm support vertical 6600-56, strong model incl. Fastening Straps 5.890-VE (2 pcs.)	visco-elastic 93.279-VE	Length 2.275-VE 560 mm	6 kg	5.750-VE	Height: 25-35 mm Width: 9-10 mm	5.40 kg