



Articulated Arms for Healthcare




FISSO
Swiss Made

Articulated Arms For Spinal Surgery

Easy, fast, accurate and reliable positioning of retractors, ring and port systems for various spine procedures

One clamping knob allows simultaneous tightening / loosening of all three joints



3D articulated arm allows pinpoint positioning of retractors

Features:

- Wide range of articulated arms (sizes, holding forces, knob colors, labeling options etc.) and peripheral components such as clamps, columns, distal holders, quick release system etc.)
- Many customer-specific designs possible
- Steam sterilizable at 134°C
- Mechanical central clamping mechanism for maintenance-free use and long lifespan
- Swiss Made





Quick release system enables fast exchange of different port sizes



Fast and stable fixation of ring systems

Our customers confirm that we are the leading company for superior quality articulated arms with central locking mechanism:



“ The FISSO articulating arms are a critical component to the success of our less invasive spine procedures. NuVasive leverages these best-in-class articulating arms to hold a variety of retractors in place securely.

The versatility, range of motion, rigidity, functionality and ease of use allows surgeons and hospitals to customize their FISSO arm set-up depending on the clinical needs of the patient which further streamline and help improve reproducibility. Without FISSO articulating arms, we would not be able to perform these procedures as routinely and reproducibly. ”

115 Hurley Road
 Unit 8D
 Oxford, CT 06478
 Telephone 860.274.5900
www.novidsurgical.com
info@novidsurgical.com

