

Patient Positioning / Arm Support System



Supine Position



Beach Chair Position



Side Position



Prone Position

Features:

- « one step set up »
- Central clamping mechanism for easy, fast and secure positioning of the arm support
- Rapid locking base allows connection to the side rails of the OP table
- Uncompromising, efficient handling and safe holding
- Space saving on the OR table, compact design
- Anatomically shaped visco-elastic (VE) – pad



Patient Arm Support Lateral



6646.03

| Order no. | Description | Pad Order no. | Articulated arm Order no. | Carrying capacity | Rail Clamp Order no. | Rail dimension | Weight |
|----------------|--|-------------------------|---------------------------|-------------------|----------------------|------------------------------------|---------|
| 6646.03 | Arm support lateral 6600-46, strong model incl. Fastening Straps 5.890-VE (2 pcs.) | visco-elastic 93.279-VE | Length 2.270-VE 460 mm | 10 kg | 5.751-VE | Height: 25-35 mm Width: 9-10 mm | 4.30 kg |

Subject to modifications

Patient Arm Support Vertical



6656.04

| Order no. | Description | Pad Order no. | Articulated arm Order no. | Carrying capacity | Rail Clamp Order no. | Rail dimension | Weight |
|----------------|---|-------------------------|---------------------------|-------------------|----------------------|------------------------------------|---------|
| 6656.04 | Arm support vertical 6600-56, strong model incl. Fastening Straps 5.890-VE (2 pcs.) | visco-elastic 93.279-VE | Length 2.275-VE 560 mm | 6 kg | 5.750-VE | Height: 25-35 mm Width: 9-10 mm | 5.40 kg |

Subject to modifications