



Articulated Arms for Healthcare

**Universal Radial Clamp for US, UK and European Surgical
Table Rails**



Easy to mount. Easy to use. Exceptional strength.



The Universal Radial Clamp

- Infinite rotation
- Fits all table rails
- Accepts 16 - 20mm accessory posts
- Steam sterilizable
- Height adjustable

Order #10-308-16 Radial Clamp

for US, UK and European rails
16mm (5/8") hole for accessory post.

Order #10-308 Radial Clamp

for US, UK and European rails
20mm hole for posts 16-20mm

Novid Surgical is the proud US Distributor for Reison AB line of surgical table accessories.