


Articulated arms 134°C Sterilizable



Article #	Type	Length mm	Bottom thread	Top thread	Knob	Weight kg
4.051-VE	XS-13	130	M6x1	M6x1	Stainless steel	0.25
4.060-VE	XS-19	190	M6x1	M6x1	Stainless steel	0.30
2.860-VE	3800-28	280	M12x1	M8x1.25	Aluminium, red	0.80
2.760-VE	3800-40	400	M12x1	M8x1.25	Aluminium, red	0.90
2.790-VE	3800-55	550	M12x1	M8x1.25	Aluminium, red	1.00
3.530-VE	4800-40	400	M12x1	M8x1.25	Aluminium, red	1.85
3.520-VE	4800-55	550	M12x1	M8x1.25	Aluminium, red	2.00
2.101-VE	5800-50	560	M12x1.75	M10x1.5	Plastic, black	3.45